Respira sicuro con

**Epidemie?
No grazie!**

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EN

Safeguard your professional image and your prestige, by protecting your health and that of you colleagues and patients. **AIRMED-BM** is the first air sterilization system.

is an air sterilizing appliance which guarantees the elimination of all microorganisms including spores. **This device can run continuously also in the presence of people, ensuring maximum safeguarding of the operators. AIRMED-BM** is simple, silent and ergonomic.

Functioning

its functioning is based on a system of forced ventilation closed loop. The air drawn into the module first passes through a filter, mechanically positioned, within the inlet entrance. Here the coarser pollutants are blocked, preventing the fouling of the germicide lamps. Then the air is forced to pass in direct contact with mercury vapor tubes that, thanks to the emission of the UV-C, carry out the maximum germicidal action. The high reflectivity shielding action of the housing concentrates the reflections of the UV-C radiation. The air is subsequently expelled through the nozzle output, thus leading to an effective microbial killing.

The system

The great advantage consists in the absolute lack of danger for the operator since there is the minimum emission of UV-C radiations from the unit. The device can effectuate a continuous air disinfection treatment (120 m³/h). It's possible to carry out a demolition of microbiologic pollution of the air in every room, during every productive phase without any damage for operator's health. **After about three hours from installation 90% of demolition of bacterial charge in the air.**

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